

<b>OrdoLog - IHE Integration Statement</b>		<b>Date</b>
		07/02/2008
<b>Vendor</b>		
<b>Product Name</b>		
<b>Version</b>		
OrdoLog	OrdoDICOM XT	2.0
<p><i>This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:</i></p>		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Scheduled Workflow	Image Manager / Image Archive Image Display Performed Procedure Step Manager	
Patient Information Reconciliation	Image Manager / Image Archive Performed Procedure Step Manager	
Access to Radiology Information	Image Manager Report Reader	
Simple Image and Numeric Report	Report Reader	
<p align="center"><b>Links to Standards Conformance Statements for the Implementation</b></p>		
<b>HL7</b> <a href="http://www.ordolog.com">www.ordolog.com</a> <b>DICOM</b> <a href="http://www.ordolog.com">www.ordolog.com</a>		
<p align="center"><b>Links to General Information on IHE</b></p>		
In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>		